

## Programmable Ramp Temp./Humid. Chamber



### ■ SPECIFICATIONS

Model	MHS
Temperature range	-20°C ~+100°C / -40°C ~+100°C / -70°C ~+100°C
Humidity range	20~95%R.H
Heating up/Pull down time	5~20°C /min
Temperature uniformity	≤±2.0°C
Interior dimensions(WxDxH)cm	80x95x70
Exterior dimensions(WxDxH)cm	136x188x225
Interior material	Stainless steel plate
Exterior material	Stainless steel plate
Volume	532L
Insulation	Rigid polyurethane foam and glass wool
Refrigeration System	Half hermetic compressor+Hermetic compressor